

#### **Contact details:**

Company: Date:

Address: Order number:

Zip code: Contact:

City: Telephone:

Country: E-mail:

# 1. Do you have a hub motor / pushing aid?

Yes go to step 9


No go to step 2

## 2. Do you want cutouts in your spokeguards?

Yes (A)

No (B)







# 3. What kind of fixation do you have?

Round mounting (A)

Flat mounting (B)





## 4. Number of fixations / cutouts

4

5

6




## 5. Dimension spokeguard - outer diameter

Outer diameter in mm: mm (max 620 mm)

For spokeguard without cutouts, see picture A.


For spokeguard with cutouts, see picture B.



### 6. Dimension spokeguard - inner diameter

Inner diameter in mm:

mm (minimum 80 mm)



#### Note:

When you remove your wheel out of your wheelchair using a quick release axle, you put your fingers through the spokes and press the axle button with your thumb.

You should have some extra space in the centre of the spokeguard to be able to press the quick release button.


# 7. Do you have a drum brake?

Yes (we recommend to choose for extra bending in step 8)


No

### 8. Do you want extra bending in your spokeguards?

Yes No









## 10. Printing

No printing, transparent spokeguards.

Design nr.: (Choose the number of the

requested design from our catalogue or our website).

Designs starting with "B-" can only be choosen for our standard dimensions without cutouts




## 11. Mouting the spokeguard on your wheel


In order to mount the spokguard on your wheel spokenipples (A) or tyraps (B) can be used. Both (A & B) are included.

If your wheel contains a hub motor, you don't have to make this choice.

These spokeguards do have standard fixation holes.







#### Choose one of these options below to determine the mounting holes

- (A) Spokeguards delivered with 3 mounting holes Ø9 mm for nipples, only to be used with spokeguard without cutouts and wheels with 36, 24, 18 or 12 spokes.
- (B) Spokeguards deliverd with 3x 2 mounting holes Ø4 mm for tyraps.
- (C) Spokeguards deliverd without any mounting holes, you will drill the mounting holes yourself.